



Chair **Caddie**TM

RETRO-FIT



It Turns...It Rolls...and LOCKS for Safety!

- Extends life of chairs by 2-4 years
- Eliminates damage to floors
- Attaches to existing chairs



The Chair CADDIE enables a spouse / caregiver to easily and safely move a seated person up-to and away-from the table...without risk of injury.

- Reduces caregiver fatigue
- Reduces workplace injuries
- Provides dignity to persons with limited mobility

NO Pushing...NO Shoving...NO Twisting!



Chair Caddie™


RETRO-FIT

1. **Attaches to existing chairs** (wood or metal)
2. Low profile design – chair is raised $\frac{3}{4}$ " above the floor
3. Fits inside the footprint of an existing chair
4. Medical grade casters:
Chair Caddie A (3" front / 2" rear)
Chair Caddie B (4" front / 3" rear)
5. Anti-tip glides
6. Front wheel tension control
7. Foot activated locking rear caster
8. Rated for 450 lbs.
9. Attaches to arm or armless chairs
10. Adjustable - Custom sizes available



Chair Caddie A

- HomeCARE
- Residential



X1 15 $\frac{3}{8}$ " - 23 $\frac{5}{8}$ "


Y 14 $\frac{3}{4}$ " - 23"

Z 16 $\frac{7}{16}$ " - 22 $\frac{9}{16}$ "

X2 10 $\frac{11}{16}$ " - 18 $\frac{15}{16}$ "

Chair Caddie B

- Senior Living / LTC
- Adult Day Care
- Commercial



X1 17 $\frac{5}{8}$ " - 24"

Y 17" - 24"

Z 18" - 24 $\frac{1}{16}$ "

X2 13" - 19 $\frac{3}{8}$ "



Watch video at www.comfortek.com/CC

Changing the *costs* culture *community* of care



www.comfortek.com | 1-888-678-2060



Watch product videos at
www.comfortek.com/better



DINING



DINING



ACTIVITY ROOM



BFL-2(1/1)



BFL-4(4/1)



BFL-4(2/2)

Butterfly-Application-Matrix-11/19

➤ *Meet Paula...*

Paula has been a resident of a care community for several years and spends most of everyday seated in a GERI Chair.

At mealtimes, Paula is transported in her GERI Chair to the dining room. Paula is positioned near a counter where care staff assist with one on one meal service, feeding her by hand as none of the furniture has been built to accommodate her chair.

Recently the community requested a Butterfly Table be set up for evaluation. Staff suggested the table be adjusted so that one of the tops was positioned suitably for Paula's use.

When Paula was brought into the dining room and moved up to her own table setting, staff were overwhelmed to see her "tear up" as for the first time she was able to feed herself with minimal assistance.

With an appropriately positioned flat surface available, staff realized there were other activities which they could do with her now that she had access to her own table.



"We had a gentleman, 6'2", sitting at a table with a lady barely 5' tall. Both, were struggling to eat on their own as the height of the table was not suitable for either of them.

With the Butterfly Table we were able to adjust the height of each table top to their specific need...further enhancing their independence."



Darell B.

Contessa



Watch product videos at
www.comfortek.com/t2-healthcare



Contessa 100



Contessa 200



Contessa 300



Contessa 100 w/ CCB
It Turns...It Rolls...and Brakes for Safety



Contessa 200 w/ CCB
It Turns...It Rolls...and Brakes for Safety



Contessa 300 w/ CCB
It Turns...It Rolls...and Brakes for Safety

Contessa-T2-Healthcare-Matrix-11/21



SENIOR LIVING / LONG TERM CARE

www.comfortek.com/t2-healthcare
1.888.678.2060